

SAFETY STOCK UP



ALL ITEMS LIMITED TO STOCK ON HAND!



3 Ply Blue Pleated Mask

- Elastic ear strap design and built in nose clip
- Not intended for surgical or certified medical use

Item #	Per Box	Price per Mask
WRBLUEPLEATMASK-10	10 Masks	\$0.40 Each



White Pleated Mask, Reusable

- Spun bond polypropylene 6 layer masks are softer, more comfortable, and more durable, than standard blue pleated paper masks
- Pliable and comfortable due to material being ultra-sonically bonded
- With metal nose clip
- Machine washable, air dry

Item #	Per Box	Price per Mask
CADMASK-6LAYER	20 Masks	\$0.87 Each



LEMYN ORGANICS™

It's What's Inside that Really Counts

Lemyn Organics™ Medical-Grade Hand Sanitizer Gel takes waterless washing to a whole new level. Our completely clean, no-rinse, residue-free formula exceeds FDA, CDC and WHO standards for effectively killing germs while remaining safe for sensitive skin. Effective. Fast-Drying. Never sticky. Deeply hydrating waterless wash. Perfect for daily use.

- Proven effective • Clean • Green • Hard-working • Hydrating**

Item #	Description	Price Each
HSLEMYNGEL-2OZ	2 oz Bottle	\$3.59 Each
HSLEMYNGEL-8OZ	8 oz Flip Top Bottle	\$6.49 Each
HSLEMYNGELP-8OZ	8 oz Pump Bottle	\$6.49 Each
HSLEMYNGEL-12OZ	12 oz Pump Bottle	\$8.49 Each

THE 99.9999% DIFFERENCE



Offer Valid: 1/4 - 1/31/2022

Please reference promo code LP1560 when placing your order with Customer Service. For online customers, enter LP1560 in the promo code field of the shopping cart.

Würth Louis and Company

(800) 422-4389 • www.wurthlac.com

OPT-OUT INFORMATION: You may request Würth Louis and Company not to send any future advertisements to a telephone facsimile machine or machines. Failure to comply with such a request in the form prescribed by law, within the shortest time prescribed by law not to exceed 30 days, is unlawful. If you wish to discontinue future fax advertisements from Würth Louis and Company, Please send your request by Email to: creditdepartment@wurthlac.com, or by Fax to: 800-944-7494 or by Calling: 800-444-0043 option 2. Your request should specify the telephone numbers or numbers of each telephone facsimile machine to which the request relates.